

Instructions for use - Please read carefully!

Lysoformin® Plus

Concentrate **Disinfecting and cleaning concentrate for surfaces, equipment and medical devices**

for professional use

Composition

100 g **solution** contains:

8.9 g Polyhexamethylene biguanide . HCl, 23.2 g Bis-(3-aminopropyl)dodecylamine,

5.0 g Didecyl Dimethyl Ammonium Chloride

Contains: (R)-p-Mentha-1,8-diene

Range of application

- **One-step** application for the **disinfection** and **cleaning** of surfaces, and equipment (incl. medical devices).
- Use according to the preparation instructions of the manufacturer of the medical devices.
- Control of nosocomial infections
- Disinfection and hygienic cleaning of all sanitary equipment.

Mode of application / dosage

Lysoformin® Plus is a liquid concentrate with which any use-dilution required can be prepared by simply adding water (at max. room temperature). Always fill container with tap water first, then add disinfectant.

Disinfection by wiping according to the VAH procedure 2002, dirty conditions		
bactericidal, levurocidal	2.00%	5 min
	1.00%	15 min
	0.75%	30 min
	0.50%	1 h
	0.25%	4 h
Tests according to the directive of the DVV/RKI		
virucidal against envelopped viruses * (incl. HIV, HBV, HCV) according to the recommendations of RKI 1/2004*	0.50%	15 min
Papova virus /Polyomavirus /SV40	0.50%	30 min
EN-Norms (low burden)		
EN 14476 (Rotavirus)	0.50%	5 min
EN 14476 (MNV)	2.00 %	2 h

*DVV 2015 (suspension test)

Lysoform – Hygiene with System

concentration	0.25 %	0.5 %	1.0%
Dosing system (for 8 l of water)			
Dosing pump for 5 l canister (1 stroke ≙ 20 ml)	1 stroke	2 strokes	4 strokes

Important: Always fill up container with water (at room temperature) and then add the disinfectant. To calculate the amount of water to be added, subtract the amount of disinfectant needed to obtain the concentration wished for from the amount of solution ready for use which has to be prepared.

Warnings / storage

- Keep disinfectants in a dry, cool but frost free and well ventilated place and out of reach of children.
- Wear suitable protective clothing, adequate protective gloves and eye protection when handling the concentrated product.
- Do not use beyond expiry date.

Properties

- Contains no aldehydes, no phenols, no phosphates; free of oxidizing compounds
- All-round efficacy, good compatibility with surfaces and materials
- Good cleaning power without addition of detergents even by heavy soil by organic substances like fat, albumin, serum
- Fresh scent
- Resists hard water
- Biodegradable

Special precautions

a) Incompatibilities:

Adding detergents to **Lysoformin® Plus** can impair or neutralize its disinfecting properties.

b) Counter-measures in case of improper use:

Should **Lysoformin® Plus** come into contact with the mucous membranes (eyes, throat, stomach etc.) get medical attention immediately.

Dilute product without delay: If it has been swallowed, drink **plenty** of water; after contact with the eyes, rinse thoroughly with **plenty** of water.

c) Miscellaneous

Objects/surfaces which come into contact with foodstuff must be rinsed with drinking water at the end of the contact/exposure time.

d) Use **Lysoformin® Plus** safely. Always read the label and product information before use.

Certified and mentioned in the relevant lists.

Standard packaging / transport units:

6 x 1000 ml bottles

2 x 5 litre canisters

article number:

1865

1869

Accessories:

Dosing pump for 5 l canister (1 stroke \cong 20 ml)

3135

Lysoform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Str. 133, D – 12247 Berlin
Tel.: +49-30-77992-0 Telefax: +49-30-77992-219
www.lysoform.de, e-mail: kontakt@lysoform.de

0318-05