

Instructions for use - Please read carefully!

Lysoformin®

Concentrate

Disinfection of surfaces and medical devices by wiping

For professional use

COMPOSITION

100 g solution contains: 6.0 g Formaldehyde, 1.8 g Glutaral, sodium alkyl ether sulphate

RANGE OF APPLICATIONS

- One-step **disinfection** and **cleaning** of surfaces and equipment in all cases of contamination by bacteria (incl. tb), viruses and skin fungus. Also suitable for the disinfection of textiles, mops and cloths.
- Control of nosocomial infections
- Disinfection and cleaning of floors and objects contaminated by skin fungus (bath-mats, grids, shoes and socks etc.)
- Laundry disinfection (tb incl.)
- Disinfection and hygienic cleaning of the whole sanitary equipment
- Disinfection and cleaning of solarium benches

MODE OF APPLICATION / DOSAGE

Lysoformin® is a liquid concentrate with which any use-dilution required can be made by simply adding water (at max. room temperature). Always fill container with tap water first and then add disinfectant. Rinse objects/surfaces that come in contact with food with drinking water after disinfection.

Disinfection by wiping accord. to the VAH-procedures of 2002 (including EN 16615), strong contamination		
bactericidal, levurocidal	1.50%	1 h
	0.75%	4 h
fungizid	4.50%	1 h
	3.00%	4 h
Tests according to the DVV/RKI directive		
Virucidal	6,00 %	1 h
Adenovirus	2.00%	1 h
Papovavirus/Polyomavirus/SV40	6.00%	1 h
EURO-Norms (slight contamination)		
EN 14476 (MNV)	2.00%	1 h

Lysoform – Hygiene with System

Dosing-system (for 8 l of water)	0.75 %	1.5 %
Dosing pump (1 stroke ≈ 30 ml)	2 strokes	4 strokes
Lysoform dosing devices	Adjust once only to obtain the desired concentration	

IMPORTANT:

The amount of water is calculated by subtracting the amount of disinfectant to be used from the total amount of solution wished for.

WARNINGS / STORAGE

- Keep disinfectants in a dry, cool but frost-free and well ventilated and dry place and out of reach of children!
- Wear suitable protective clothing, adequate protective gloves and eye protection when handling the concentrated product.
- Do not use beyond expiry date.
- Shelf-life of the ready-for-use solutions* without residue of organic material * and when stored properly (in a cool and dry place) in clean and tightly closed containers: 4 weeks.

* The period during which the proven efficacy of the disinfectant solution is maintained if, as recommended by the manufacturer, it has been used at the recommended concentration and thereafter covered and left unused for the period specified in the manufacturer's instructions

PROPERTIES

- Bactericidal (incl. tb), fungicidal, virucidal
- Good compatibility with materials
- Freshly scented
- Biodegradable

SPECIAL PRECAUTIONS

a) Incompatibilities

Adding detergents to **Lysoformin**[®] can impair or neutralize the disinfecting properties.

b) Counter-measures in case of improper use

Should **Lysoformin**[®] come into contact with the mucous membranes (eyes, throat, stomach, etc.) get medical attention immediately. Dilute product without delay: if it has been swallowed, drink **plenty** of water; after contact with the eyes, rinse **thoroughly** with plenty of water.

c) Use **Lysoformin**[®] safely. Always read the label and product information before use.

Certified and mentioned in the relevant lists.

STANDARD PACKAGING / TRANSPORT UNITS

Article numbers

12 x 1 litre bottle	3635
2 x 5 litre canister	3695
1 x 20 litre canister	5270

Accessories

Dosing pump for 5 litre canister (1 pump stroke = 30 ml)	3136
---	------

Dosing system

Lysoform dosing device DH 9 I	9174
-------------------------------	------

Lysoform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Str. 133, D – 12247 Berlin
Tel.: +49-30-77992-0 Telefax: +49-30-77992-219
www.lysoform.de, e-mail: kontakt@lysoform.de

0218-11