

Instructions for use - Please read carefully!

# Fugaten®-Spray

(Lysoform®-Spray)

Ready for use

Disinfection-Spray based on alcohol for the rapid disinfection of surfaces, objects and medical devices

For professional use

## Composition

100 g solution contains: 70 g Ethanol (denatured)

## Range of application

For rapid disinfection of

- Small surfaces, objects and medical devices (e.g. operating tables, examination chairs etc.)
- Surfaces difficult of access
- All alcohol resistant materials
- Product exists also without perfume, particularly adequate for use in the food industry

Disinfection by wiping according to the VAH procedure 2015 mechanical effect, clean conditions		
bactericidal, levurocidal	100%	60 s
Tests according to the directive of DVV/RKI		
virucidal against enveloped viruses (incl. HIV, HBV, HCV) according to the recommendations of RKI 1/2004	100%	30 s
rota virus	100%	30 s
adeno virus	100%	30 s
papova virus/polyoma virus/SV 40	100%	15 mn
EN-Norms		
EN 14348 (tuberculicidal, mycobactericidal), dirty conditions	80%	30 s
EN 14476 (polio virus), clean conditions	100%	5 mn
EN 14476 (MNV), clean conditions	100%	30 s
virucidal	100%	15 mn

## Warnings / storage

- **Fugaten®-Spray** is easily inflammable!
- Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. Do not smoke.  
Max. quantity to be sprayed: 50 ml/ m<sup>2</sup> of the surface of the room
- Keep disinfectant in a dry, cool but frost-free and well-ventilated area, away from ignition sources and out of reach of children.
- Do not use beyond expiry date.

## Properties

- Bactericidal (incl. tb), fungicidal, virucidal
- Rapid action and drying after spraying on surfaces and objects
- Freshly scented; deodorizes the surrounding air
- Good wetting power
- No adverse effect on most materials
- Ready to use product in refillable, polyethylene bottles with spray-pump

### **Special precautions**

a) Counter-measures in case of improper use

**Fugaten®-Spray** is irritating in contact with the mucous membranes (eyes, throat, stomach, bowels etc.). In case of such irritations dilute product immediately. If it has been taken **orally**, make patient **drink plenty** of water; if it has been in contact with other mucous membranes, rinse thoroughly with **plenty** of water. In each case, get medical attention.

b) Miscellaneous

Do not spray on materials which are not compatible with alcohol (coated, pickled surfaces, Plexiglas etc.).

c) Use disinfectant safely. Always read the label and product information before use.

Certified and mentioned in the relevant lists.

### **Standard packaging/transport units:**

10 x 500 ml bottle with spray pump

12 x 1000 ml bottle with spray pump

12 x 1000 ml bottle with spray pump (without perfume)

2 x 5 litre canister

2 x 5 litre canister (without perfume)

**Lysoform Dr. Hans Rosemann GmbH**  
Kaiser-Wilhelm-Str. 133, D – 12247 Berlin  
Tel.: +49-30-77992-0    Telefax: +49-30-77992-219  
[www.lysoform.de](http://www.lysoform.de), e-mail: kontakt@lysoform.de

**1118-15**